**GBRUIKSAANWIJZING  
MODE D'EMPLOI  
USERS MANUAL  
GEBRAUCHSANWEISUNG  
MANUAL DE INSTRUCCIONES  
MANUAL DE INSTRUÇÕES**

**BC-02  
TWIN THERMOMETER SET**  
**Alecto®**

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**GB General:**  
The Alecto thermometer set allows you to accurately and conveniently measure the temperature of the human body. The set contains the following items:

**Model 1**, for adults and grown up children: for rectal use (in the anus), for oral use (in the mouth) or for under the arm

**Model 2**, special for baby's: for oral use (in the mouth) (measure after 6-8 min)

Please read this guide carefully before using one of the thermometers and keep it for future reference.

**Important safeguards:**

- Never attach any string or cord to soother-thermometer. It could strangle the child.
- Before each use check the nipple of the soothers for tears or cracks and pull it to ensure that it is firmly attached to the shield. Do not use it if a defect is found or if it is loose.
- Do not use the soother thermometer as a regular soother
- Clean the thermometers every time before and after use.
- Store the thermometers out of reach of children.

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**Memory:**  
As soon as you switch-on the thermometer, the last measured temperature comes on the LCD, together with a M. Attention: this readout is on the LCD only momentarily and only once per measurement.

**Cleaning:**  
To clean the thermometers after each use, wipe it with a damp cloth only or wash the unit with mild soap and warm water. The thermo-meters are water resistant only, not waterproof.

Never sterilize nor boil the thermometers, do not wash in dishwasher, do not put in microwave or microwave sterilizer and never use any chemical solvents.

Tip: use special plastic hygiene covers to cover the tip of the thermometer, these are available from most pharmacies.

**Battery:**  
The thermometers have long-life batteries which should last for about 12000 min.

**Replacing battery (model-1):**  
1 Remove the battery cover.  
2 Use an insulated probe to remove the battery

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3 Insert the new battery (LR41) with the positive side (+) facing upwards.

4 Secure the cover back in place.

Battery: 1.5V, LR41 alkaline or equivalent.

**Model 2:**  
The battery cannot be replaced. Dispose of the thermometer in an environmentally friendly manner by hand it in at your supplier.

**Battery save mode:**  
In case you do not manually switch-off the thermometer, it switches off automatically after 10 minutes.

**Caring:**

- Don't drop thermometers or expose them to shocks or high temperatures.
- Only use these thermometers for the measurement of human body temperatures.
- Do not store in direct sunlight or near a heat source.
- Do not bend or disassemble the thermometers.

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